

## Parker Tub Cushioning and Footrest Options

### 1. Cushioning options to prevent patients from slipping forward:



**Moderate Anti-slip Needs**

Parker tub headrest cushion positioned at front of seat pan



**Significant Anti-slip Needs**

Gray headrest pillow positioned at front of seat pan



**Significant Anti-slip Needs**

Velcro thigh strap attached around the side grab bars and above patient's legs with headrest cushion beneath

### 2. Cushioning options for patients' needs:



Cushioning option for **patients who require increased neck support**



Cushioning option for **patients who require increased lower back support**



Cushioning option for **patients who require increased support on one side (weakness on one side)**

## 2. Cushioning options for patients' needs (continued):



Cushioning option for patients who require increased support on both sides (slimmer bodies)



Cushioning options for patients with coccygeal ulcers; also provides increased access to the person's perineal area



Cushioning options for patients with a frail stature and/or a bony spine

## 3. Footrest options:



Footrest option for staff to give foot care to patient in tub.



Can use footrest at the top of the tub as patients may feel more secure with this option. Soap and washcloth can be placed on footrest as well.