

Ver 1.0 Sep 14, 2009

COLLECTION OF SALIVA FOR CORTISOL

Follow these instructions exactly. The accuracy of test results is dependent on the quality of the specimen collected.

- 1) **DO NOT** eat or drink 30 minutes prior to collection of saliva.
- 2) **DO NOT** brush teeth or rinse mouth 30 minutes prior to collection of saliva.
- 3) Label collection container with patient's name, date of birth, collection date and the **time using the 24-hour clock e.g., 23:00 or 24:00.**

Collect specimen between 23:00 and 24:00 unless otherwise directed (for Cushing's syndrome, midnight saliva is the specimen to collect).

- 4) At designated time remove the TOP cap of the plastic container to expose the swab.
- 5) Tip the container in order to place the swab in mouth. **DO NOT TOUCH THE SWAB WITH FINGERS.**
- 6) Hold the swab in the mouth for 2 full minutes or longer **to ensure sufficient volume of sample collected (minimum 0.6 mL, 2 mL preferred)**. Roll it around with the tongue, but **DO NOT CHEW OR SUCK ON THE SWAB.**
- 7) Place the swab back into the plastic container in the same location it came from. **DO NOT TOUCH THE SWAB WITH FINGERS.** Replace the cap.
- 8) Store the sample in the refrigerator until you can deliver it to the lab. Ensure you have recorded the time, **23:00 or 24:00** and date on the sample.
- 9) **If you drop off your salivette tube at lab other than VGH, please instruct the lab to forward the collection container to the address below.**

SEND TO:

Vancouver General Hospital

Laboratory Reception

910 West 10th Avenue Room 1300 JPPN

Vancouver BC V5Z 4E3

Canada

Tel: 1-877-747-2522 Canada

1-604-875-4111 (54811) Outside Canada

Ver 1.0 Sep 14, 2009

