#### 828468 Feb 7-17

Healthlink BC

Call 8-1-1

www.healthlinkbc.ca

# How to use a Turbuhaler®

#### Patient Information

### Usage

- You do not need to shake the Turbuhaler<sup>®</sup> before using. If you do happen to shake it you will hear a sound. This sound is drying agent, not the medication.
- If you breathe into the Turbuhaler<sup>®</sup> or drop it after it is loaded, you will need to load another dose. Follow the steps from the beginning. You won't get a double dose.

#### Cleaning

- Wipe the mouthpiece with a clean, dry, cloth or tissue.
- Never use water.

#### Storage

- Keep the cap on when not in use.
- Keep it in a dry place at room temperature.
- Keep away from direct sunlight or heat.

Nurse

Dietitian

Pharmacist

Hearing Impaired

24 hours a day

9 am-5 pm

5 pm-9 am

Call 7-1-1

Daily

Daily

Mon-Fri

#### When to get a new Turbuhaler®

- Turbuhaler<sup>®</sup> has a window that shows you when the inhaler is empty and there are 2 types. Please see the package insert to find out when your inhaler is empty.
- Check the expiry date on your inhaler.
- You should have a new inhaler ready to use when the old one is empty or expired.

This document is for informational purposes only. For detailed instructions on use of the Turbuhaler<sup>®</sup> and information about your medication please see the package insert or ask your pharmacist.

more information on other side  $\rightarrow$ 

Call 8-1-1 to speak with a nurse,


ask a dietitian about nutrition, or a

pharmacist about your medication.

Interior Health Every person matters

www.interiorhealth.ca

## How to use a Turbuhaler $^{\ensuremath{\mathbb{R}}}$ (continued)



**I. Hold** the Turbuhaler<sup>®</sup> with the coloured grip on the bottom. **Twist** the cap off.



**2.** Turn the grip one way and back. You will hear a click.



**3.** Breathe out, away from the device.



**4.** Put the Turbuhaler<sup>®</sup> between your teeth and seal your lips around it. **Breathe in very fast,** until your lungs are full.



5. Remove the Turbuhaler<sup>®</sup> and hold your breath for 10 seconds, if you can. Breathe out. If you need to take another puff repeat steps 2 to 5.



**6.** Twist the cap on.



7. Rinse, gargle and spit out.