How to use a Diskus®



Patient Information www.interiorhealth.ca

Usage

- Do not shake the Diskus® before using.
- If you drop or breathe into the Diskus® after loading it (Step I and 2), you will need to load another dose. Close the Diskus® and follow the steps from the beginning. Don't worry you won't get a double dose.

Cleaning

- Wipe the mouthpiece with a clean, dry, cloth or tissue after you take your dose.
- Never use water.

Storage

- Keep it closed.
- Keep it in a dry place at room temperature.
- Keep away from direct sunlight or heat.

When to get a new Diskus®

- Check the dose counter. When the window of the dose counter is red and you see the number '0' your Diskus[®] is empty.
- Check the expiry date on the Diskus[®].
- Start a new inhaler when the old one is empty or expired.

This document is for informational purposes only. For detailed instructions on use of the Diskus® and information about your medication please see the package insert or ask your pharmacist.

more information on other side >


Healthlink BC	Nurse	24 hours a day	Daily
Call 8-1-1	Dietitian	9 am-5 pm	Mon-Fri
	Pharmacist	5 pm-9 am	Daily
www.healthlinkbc.ca	Hearing Impaired	Call 7-1-1	

Call 8-1-1 to speak with a nurse, ask a dietitian about nutrition, or a pharmacist about your medication.

How to use a Diskus® (continued)



• Hold the Diskus® with your thumb in the groove. Push Open. Hear a click.



2. Slide lever. Hear a click.



3. Breathe out away from the Diskus[®] to the end of a normal breath.



4. Put the Diskus® between your teeth and seal your lips around it. Breathe in as fast as you can until your lungs are full.



5. Remove the Diskus® and hold your breath for 10 seconds, if you can. Breathe out.



6. Close.



7. Rinse, gargle and spit out.